

## Pediatric Three-Lumen CVC for Femoral Vein Insertion

**CVC SET: 3L 5.5 FR X 30 CM**

SKU / Article #: CS-14553

### Highlights

CE Mark

### Considerations

Not made with natural rubber latex.

**Sales Quantity/Case: 5**





*Representative photo(s): actual product specifications may vary from that shown in the photo*

#### Pediatric Three-Lumen CVC for Femoral Vein Insertion

SKU / Article #: CS-14553

Quantity	Component Description
1	Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 30 cm
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
1	Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
1	Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
3	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.