

## Pediatric Two-Lumen CVC

### CVC SET: 2-LUMEN 4 FR X 8 CM

SKU / Article #: CS-15402-E

### Highlights

CE Mark

### Considerations





Not made with natural rubber latex.

**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Pediatric Two-Lumen CVC SKU / Article #: CS-15402-E	
Quantity	Component Description
1	Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 8 cm
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
1	Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
1	Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.