

Pediatric Three-Lumen CVC

CVC SET: 3-LUMEN 5.5 FR X 8 CM

SKU / Article #: CS-15553-E

Highlights

CE Mark

Considerations





Not made with natural rubber latex.

Sales Quantity/Case: 5

Pediatric Three-Lumen CVC SKU / Article #: CS-15553-E	
Quantity	Component Description
1	Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 8 cm
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
1	Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
1	Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
3	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp

Representative photo(s): actual product specifications may vary from that shown in the photo

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.