

Four-Lumen CVC

CVC SET: 4-LUMEN 8.5 FR X 20 CM

SKU / Article #: CS-15854-E

Highlights

CE Mark

Considerations



Not made with natural rubber latex.

Sales Quantity/Case: 5

Four-Lumen CVC SKU / Article #: CS-15854-E	
Quantity	Component Description
1	Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 20 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Pressure Transduction Probe
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
1	Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
4	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp

Representative photo(s): actual product specifications may vary from that shown in the photo

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.