

## Arrowgard Blue® Two-Lumen Hemodialysis Curved Catheter for High Volume Infusions

HEMODIALYSIS SET: 2-LUMEN 14 FR X 15  
CM

SKU / Article #: CS-22142-CF

### Highlights

---

Antimicrobial: Arrowgard Blue\*  
CE Mark

### Considerations

---

\*Chlorhexidine  
Not made with natural rubber latex.

**Sales Quantity/Case: 5**

---




*Representative photo(s): actual product specifications may vary from that shown in the photo*

# Arrow+ard Blue® Two-Lumen Hemodialysis Curved Catheter for High Volume Infusions

SKU / Article #: CS-22142-CF

Quantity	Component Description
1	Two-Lumen Catheter: Curved: 14 Fr. (4.72 mm OD) x 15 cm Radiopaque Polyurethane with Blue FlexTip®, Arrow+ard® Antimicrobial Surface Treatment, Extension Line Clamps
1	Spring-Wire Guide, Marked: .038" (0.97 mm) dia. x 27-1/2" (70 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 14 Fr. (4.7 mm) x 14 cm
1	Tissue Dilator 16 Fr. (5.3 mm) x 4 1/2" (11.4 cm)
2	Dust Cap: Non-Vented
1	Suture Wing
1	Drape: 24" x 36" (61 cm x 91 cm)
1	Scalpel: #11

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.