

Pediatric Arrowg+ard Blue[®] Two-Lumen CVC

CVC SET: 2-LUMEN 5 FR X 5 CM

SKU / Article #: CS-22502

Highlights

Antimicrobial: Arrowg+ard Blue*
CE Mark

Considerations




*Chlorhexidine
Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Pediatric Arrowg+ard Blue [®] Two-Lumen CVC	
SKU / Article #: CS-22502	
Quantity	Component Description
1	Arrowg+ard Blue [®] Two-Lumen Catheter: 5 Fr. (1.8 mm OD) x 5 cm
1	Spring-Wire Guide: .021" (0.53 mm) dia. x 17-3/4" (45 cm) PTFE-Coated (Straight Soft Tip on One End - "J" Tip on Other)
1	Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
1	Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle with 5 mL Luer-Slip Syringe
1	Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
1	Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
2	Dust Cap: Non-Vented

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.