

Representative photo(s): actual product specifications may vary from that shown in the photo

Arrowg+ard Blue® Two-Lumen CVC

CVC SET: 2-LUMEN 8FR X 16CM

SKU / Article #: CS-22802

Highlights

Antimicrobial: Arrowg+ard Blue*
CE Mark

Considerations




*Chlorhexidine
Not made with natural rubber latex.

Sales Quantity/Case: 5

Arrowg+ard Blue® Two-Lumen CVC
SKU / Article #: CS-22802

Quantity	Component Description
1	Arrowg+ard Blue® Two-Lumen Catheter: 8 Fr. (2.8 mm OD) x 16 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.