



## Arrowgard Blue® Four-Lumen CVC

CVC SET: 4-LUMEN 8.5 FR X 16 CM

SKU / Article #: CS-22854-E

### Highlights

Antimicrobial: Arrowgard Blue\*  
CE Mark

### Considerations

\*Chlorhexidine  
Not made with natural rubber latex.

**Sales Quantity/Case: 5**

Arrowgard Blue® Four-Lumen CVC SKU / Article #: CS-22854-E	
Quantity	Component Description
1	Arrowgard Blue® Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 16 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Pressure Transduction Probe
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
1	Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
4	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp

Representative photo(s): actual product specifications may vary from that shown in the photo

