

## Arrowgard Blue® Three-Lumen CVC

CVC SET: 3-LUMEN 12 FR X 20 CM

SKU / Article #: CS-25123-F

### Highlights

Antimicrobial: Arrowgard Blue\*  
CE Mark

### Considerations





\*Chlorhexidine  
Not made with natural rubber latex.

**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Arrowgard Blue® Three-Lumen CVC SKU / Article #: CS-25123-F	
Quantity	Component Description
1	Arrowgard Blue® Three-Lumen Catheter: 12 Fr. (4.1 mm OD) x 20 cm
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Pressure Transduction Probe
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
1	Drape: 24" x 36" (61 cm x 91 cm)
1	Scalpel: #11

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.