

Pediatric Arrowg+ard Blue® Three-Lumen CVC

CVC SET: 3-LUMEN 5.5FR X 13CM

SKU / Article #: CS-26553

Highlights

Antimicrobial: Arrowg+ard Blue*
CE Mark

Considerations




*Chlorhexidine
Not made with natural rubber latex.

Sales Quantity/Case: 5

Pediatric Arrowg+ard Blue® Three-Lumen CVC	
SKU / Article #: CS-26553	
Quantity	Component Description
1	Arrowg+ard Blue® Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 13 cm
1	Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
1	Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
3	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp

Representative photo(s): actual product specifications may vary from that shown in the photo

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.