

## Arrowgard Blue Plus® Two-Lumen CVC

### CVC SET: 2-LUMEN 8 FR X 16 CM

SKU / Article #: CS-42802

#### Highlights

Antimicrobial: Arrowgard Blue Plus\*  
CE Mark

#### Considerations

\*Chlorhexidine  
Not made with natural rubber latex.




**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Arrowgard Blue Plus® Two-Lumen CVC	
SKU / Article #: CS-42802	
Quantity	Component Description
1	Arrowgard Blue Plus® <sup>1</sup> Two-Lumen Catheter: 8 Fr. (2.8 mm OD) x 16 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp

<sup>1</sup> Licensed under US Patent No. 7,329,412.

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.