



Representative photo(s): actual product specifications may vary from that shown in the photo

You-Bend™ Arrowg+ard Blue® Two-Lumen Hemodialysis Catheter for High Volume Infusions

HEMODIALYSIS SET: 2-LUMEN 12 FR X 20 CM

SKU / Article #: CU-25122-F

Highlights

Antimicrobial: Arrowg+ard Blue\*  
CE Mark


Considerations

\*Chlorhexidine  
Not made with natural rubber latex.  
‡MR Conditional

Sales Quantity/Case: 5

| You-Bend™ Arrowg+ard Blue® Two-Lumen Hemodialysis Catheter for High Volume Infusions<br>SKU / Article #: CU-25122-F |   |
|---|---|
| Quantity  | Component Description   |
| 1   | Two-Lumen Catheter‡: 12 Fr. (4.06 mm OD) x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowg+ard® Antimicrobial Surface Treatment, Extension Line Clamps |
| 1   | Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer                         |
| 1   | Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle   |
| 1   | Pressure Transduction Probe   |
| 1   | Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe   |
| 1   | Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm   |
| 2   | Dust Cap: Non-Vented  |
| 1   | Drape: 24" x 36" (61 cm x 91 cm)  |
| 1   | Scalpel: #11  |

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.