

Representative photo(s): actual product specifications may vary from that shown in the photo

You-Bend™ Two-Lumen Hemodialysis Catheter for High Volume Infusions

HEMODIALYSIS SET: 2-LUMEN 12 FR X 13 CM

SKU / Article #: CV-13122-UF

Highlights

CE Mark





Considerations

Not made with natural rubber latex.
#MR Conditional

Sales Quantity/Case: 5

You-Bend™ Two-Lumen Hemodialysis Catheter for High Volume Infusions SKU / Article #: CV-13122-UF	
Quantity	Component Description
1	Two-Lumen Catheter‡: 12 Fr. (4.06 mm OD) x 13 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
2	Dust Cap: Non-Vented

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.