

## Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

**CVC SET: 3-LUMEN 7 FR X 20 CM**

SKU / Article #: CV-15703-M

### Considerations



Not made with natural rubber latex.

**Sales Quantity/Case: 10**

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter	
SKU / Article #: CV-15703-M	
Quantity	Component Description
1	Multi-Lumen Indwelling Catheter: 7 Fr. (2.5 mm OD) x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 60 cm (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 18 Ga. x 6.35 to 7.20 cm XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator
3	Dust Cap: Non-Vented
1	Fastener: Catheter Clamp
1	Clamp: Catheter

*Representative photo(s): actual product specifications may vary from that shown in the photo*

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.