

## Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

**CVC SET: 2-LUMEN 7 FR X 20 CM**

SKU / Article #: CV-17702-M

### Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 10




*Representative photo(s): actual product specifications may vary from that shown in the photo*

### Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

SKU / Article #: CV-17702-M

Quantity	Component Description
1	Two-Lumen Indwelling Catheter: 7 Fr. (2.5 mm OD) x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 60 cm (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 18 Ga. x 6.35 to 7.20 cm XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator
2	Dust Cap: Non-Vented
1	Fastener: Catheter Clamp
1	Clamp: Catheter

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.