

## Arrowgard Blue Plus® Two-Lumen CVC

### CVC SET: 2-LUMEN 7 FR X 20 CM

SKU / Article #: CV-47702-E

#### Highlights

Antimicrobial: Arrowgard Blue Plus\*  
CE Mark

#### Considerations




\*Chlorhexidine  
Not made with natural rubber latex.

**Sales Quantity/Case: 10**

Arrowgard Blue Plus® Two-Lumen CVC	
SKU / Article #: CV-47702-E	
Quantity	Component Description
1	Arrowgard Blue Plus® <sup>1</sup> Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp

<sup>1</sup> Licensed under US Patent No. 7,329,412.

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.