



## Five-Lumen CVC

### CVC SET: 5-LUMEN 9.5FR X 16CM

SKU / Article #: DE-12955-S

### Highlights

Sharps Safety  
CE Mark

### Considerations

Not made with natural rubber latex.

**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Five-Lumen CVC	
SKU / Article #: DE-12955-S	
Quantity	Component Description
1	Five-Lumen Catheter: 9.5 Fr. (3.3 mm OD) x 16 cm
1	Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
5	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	Safety Scalpel: #11

