

Representative photo(s): actual product specifications may vary from that shown in the photo

Two-Lumen CVC

CVC SET: 2-LUMEN 7FR X 20CM

SKU / Article #: DE-17702-S

Highlights

Sharps Safety  
CE Mark


Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

Two-Lumen CVC SKU / Article #: DE-17702-S	
Quantity	Component Description
1	Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm
1	Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	Safety Scalpel: #11

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.