

## Arrowg+ard Blue® Three-Lumen CVC

CVC SET: 3-LUMEN 8.5FR X 20CM

SKU / Article #: DE-25853-S

### Highlights

Antimicrobial: Arrowg+ard Blue\*  
Sharps Safety  
CE Mark

### Considerations

\*Chlorhexidine  
Not made with natural rubber latex.

**Sales Quantity/Case: 5**






*Representative photo(s): actual product specifications may vary from that shown in the photo*

# Arrowgard Blue® Three-Lumen CVC

SKU / Article #: DE-25853-S

Quantity	Component Description
1	Arrowgard Blue® Three-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 20 cm
1	Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
3	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	Safety Scalpel: #11

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.