

Arrowgard Blue® Four-Lumen CVC

CVC SET: 4-LUMEN 8.5FR X 20CM

SKU / Article #: DE-25854-S

Highlights

Antimicrobial: Arrowgard Blue*
 Sharps Safety
 CE Mark

Considerations




*Chlorhexidine
 Not made with natural rubber latex.

Sales Quantity/Case: 5

Arrowgard Blue® Four-Lumen CVC SKU / Article #: DE-25854-S	
Quantity	Component Description
1	Arrowgard Blue® Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 20 cm
1	Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
4	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	Safety Scalpel: #11

Representative photo(s): actual product specifications may vary from that shown in the photo

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.