

Pressure Injectable Three-Lumen PICC pre-loaded with Arrow® NaviCurve™ Stylet

PI PICC NaviCurve: 3L 6FR x 55CM w BP

SKU / Article #: DLX-35563-CURVB

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.
Includes Lidocaine
##MR Safe

Sales Quantity/Case: 3

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Three-Lumen PICC pre-loaded with Arrow® NaviCurve™ Stylet SKU / Article #: DLX-35563-CURVB

Quantity	Component Description
1	Three Lumen TaperFree® Catheter###: 6 Fr. (2.14 mm OD) x 55 cm, Pressure Injectable, Arrow® NaviCurve™ Stylet and T-Port Connector
1	GlideThru™ Peel-Away Sheath: 6 Fr. x 2-3/4" (7 cm) Radiopaque over 6 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Injection Needle: SafetyGlide™1 25 Ga. x 5/8" (1.60 cm) RW
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 5 mL Luer-Lock
4	10 mL Luer-Lock Pre-Filled Saline Syringe
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep

Pressure Injectable Three-Lumen PICC pre-loaded with Arrow NaviCurve Stylet
 SKU / Article #: DLX-35563-CURVB


Quantity	Component Description
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep ^{®2} One-Step Solution
2	3 mL Applicator 2% CHG and 70% IPA Chloraprep ^{®2} One-Step Solution with Hi-Lite Orange [™] Tint
3	Dust Cap: Non-Vented
1	SharpsAway [®] II Locking Disposal Cup
1	SharpsAway [®] Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape [™] with two 13 cm (5") x 15 cm (6") fenestrations
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Patient ID Card
1	Patient Information Booklet
1	Medication Label: 1% Lidocaine
2	Paper Tape Measure
1	Tourniquet
1	Foam Electrodes (3 per pack)
10	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Dressing: BIOPATCH ^{®3}
1	Transducer Cover with Gel and Rubber Bands
1	Remote Cover
1	Surgical Apparel: Mask with Eye Shield
1	Dressing: STATLOCK ^{®4} PICC Plus Catheter Stabilization Device
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask

Pressure Injectable Three-Lumen PICC pre-loaded with Arrow NaviCurve Stylet
 SKU / Article #: DLX-35563-CURVB

Quantity	Component Description
1	Tape: Steri-Strip®5
1	T-Piece Cover

- ¹ A trademark of Becton, Dickinson and Company.
- ² A registered trademark of CareFusion or one of its subsidiaries.
- ³ A registered trademark of Johnson & Johnson Corporation.
- ⁴ A registered trademark of C. R. Bard, Inc.
- ⁵ A registered trademark of 3M Company.

Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.