

Pressure Injectable Arrowg+ard Blue Advance® One-Lumen PICC pre-loaded with Arrow® VPS Precision™ Stylet

AGBA PICC G4 Stylet: 1L 4.5FR x 55CM BP

SKU / Article #: DLX-45541-VPSB

**Highlights**

- Antimicrobial: Arrowg+ard Blue Advance\*
- Antithrombogenic: Arrowg+ard Blue Advance\*
- Sharps Safety
- Pre-loaded with VPS Stylet

**Considerations**

- \*Chlorhexidine
- Not made with natural rubber latex.
- Includes Lidocaine
- ###MR Safe

**Sales Quantity/Case: 3**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Pressure Injectable Arrowg+ard Blue Advance® One-Lumen PICC pre-loaded with Arrow® VPS Precision™ Stylet

SKU / Article #: DLX-45541-VPSB

| Quantity | Component Description  |
|----------|--|
| 1        | One Lumen TaperFree® Catheter### with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection <sup>1</sup> : 4.5 Fr. (1.58 mm OD) x 55 cm, Pressure Injectable, Contamination Guard, Arrow® VPS Precision™ Stylet and Tuohy-Borst Adapter |
| 1        | GlideThru™ Peel-Away Sheath: 4.5 Fr. x 2-3/4" (7 cm) Radiopaque over 4.5 Fr. Dilator   |
| 1        | Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer  |
| 1        | Stiffening Stylet (0.36mm dia. x 75cm with straight distal tip)  |
| 1        | Injection Needle: SafetyGlide™ <sup>2</sup> 25 Ga. x 5/8" (1.60 cm) RW   |
| 1        | Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW  |
| 1        | Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)   |

Pressure Injectable Arrowgard Blue Advance One-Lumen PICC pre-loaded with Arrow VPS Precision Stylet

SKU / Article #: DLX-45541-VPSB


| Quantity | Component Description   |
|----------|---|
| 1        | Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW  |
| 1        | Syringe: 5 mL Luer-Lock   |
| 4        | 10 mL Luer-Lock Pre-Filled Saline Syringe   |
| 1        | 5 mL 1% Lidocaine HCl Solution and Alcohol Prep   |
| 1        | 3 mL Applicator 2% CHG and 70% IPA Chloraprep <sup>®3</sup> One-Step Solution                                       |
| 2        | 3 mL Applicator 2% CHG and 70% IPA Chloraprep <sup>®3</sup> One-Step Solution with Hi-Lite Orange <sup>™</sup> Tint |
| 1        | Dust Cap: Non-Vented  |
| 1        | SharpsAway <sup>®</sup> II Locking Disposal Cup   |
| 1        | SharpsAway <sup>®</sup> Disposal Cup  |
| 1        | Catheter Trimmer  |
| 1        | Towel: 24" x 36" (61 cm x 91 cm)  |
| 1        | Maximal Barrier Drape <sup>™</sup> with two 13 cm (5") x 15 cm (6") fenestrations                                   |
| 1        | Safety Scalpel: #11   |
| 1        | Sterile Procedure Sign  |
| 1        | Patient ID Card   |
| 1        | Patient Information Booklet   |
| 1        | Medication Label: 1% Lidocaine  |
| 2        | Paper Tape Measure  |
| 1        | Tourniquet  |
| 1        | Foam Electrodes (3 per pack)  |
| 10       | Gauze Pad: 2" x 2" (5 cm x 5 cm)  |
| 10       | Gauze Pad: 4" x 4" (10 cm x 10 cm)  |
| 1        | Surgical Apparel: Impervious Gown   |
| 1        | Dressing: I.V. Advanced 8.5 cm x 11.5 cm  |
| 1        | Dressing: BIOPATCH <sup>®4</sup>  |
| 1        | Transducer Cover with Gel and Rubber Bands  |
| 1        | Surgical Apparel: Mask with Eye Shield  |

Pressure Injectable Arrowgard Blue Advance One-Lumen PICC pre-loaded with Arrow VPS Precision Stylet  
 SKU / Article #: DLX-45541-VPSB

| Quantity | Component Description  |
|----------|--|
| 1        | Dressing: STATLOCK <sup>®5</sup> PICC Plus Catheter Stabilization Device |
| 1        | Surgical Apparel: Bouffant Cap   |
| 1        | Surgical Apparel: Mask   |
| 1        | Tape: Steri-Strip <sup>®6</sup>  |

- <sup>1</sup> Licensed under US Patent No. 7,329,412 and 11,219,706.
- <sup>2</sup> A trademark of Becton, Dickinson and Company.
- <sup>3</sup> A registered trademark of CareFusion or one of its subsidiaries.
- <sup>4</sup> A registered trademark of Johnson & Johnson Corporation.
- <sup>5</sup> A registered trademark of C. R. Bard, Inc.
- <sup>6</sup> A registered trademark of 3M Company.

## Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.