



Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen PICC pre-loaded with Arrow® NaviCurve™ Stylet

AGBA PICC NaviCurve: 2L 5.5FR x 55CM TCG

SKU / Article #: DLX-45552-CURVC

Highlights

- Antimicrobial: Arrowg+ard Blue Advance*
- Antithrombogenic: Arrowg+ard Blue Advance*
- Sharps Safety

Considerations

- *Chlorhexidine
- Not made with natural rubber latex.
- Includes Lidocaine
- ###MR Safe

Sales Quantity/Case: 3

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen PICC pre-loaded with Arrow® NaviCurve™ Stylet	
SKU / Article #: DLX-45552-CURVC	
Quantity	Component Description
1	Two Lumen TaperFree® Catheter### with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection ¹ : 5.5 Fr. (1.88 mm OD) x 55 cm, Pressure Injectable, Arrow® NaviCurve™ Stylet and T-Port Connector
1	GlideThru™ Peel-Away Sheath: 5.5 Fr. x 2-3/4" (7 cm) Radiopaque over 5.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Injection Needle: SafetyGlide™ ² 25 Ga. x 5/8" (1.60 cm) RW
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW

Pressure Injectable Arrowgard Blue Advance Two-Lumen PICC pre-loaded with Arrow NaviCurve Stylet
 SKU / Article #: DLX-45552-CURVC

Quantity	Component Description
1	Syringe: 5 mL Luer-Lock
4	10 mL Luer-Lock Pre-Filled Saline Syringe
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep ^{®3} One-Step Solution
2	3 mL Applicator 2% CHG and 70% IPA Chloraprep ^{®3} One-Step Solution with Hi-Lite Orange [™] Tint
2	Dust Cap: Non-Vented
1	SharpsAway [®] II Locking Disposal Cup
1	SharpsAway [®] Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape [™] with two 13 cm (5") x 15 cm (6") fenestrations
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Patient ID Card
1	Patient Information Booklet
1	Medication Label: 1% Lidocaine
2	Paper Tape Measure
1	Tourniquet
1	Foam Electrodes (3 per pack)
10	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: Tegaderm ^{®4} CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm)
1	Transducer Cover with Gel and Rubber Bands
1	Remote Cover
1	Surgical Apparel: Mask with Eye Shield
1	Dressing: STATLOCK ^{®5} PICC Plus Catheter Stabilization Device

Pressure Injectable Arrowgard Blue Advance Two-Lumen PICC pre-loaded with Arrow NaviCurve Stylet
 SKU / Article #: DLX-45552-CURVC

Quantity	Component Description
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask
1	Tape: Steri-Strip®4
1	T-Piece Cover

¹ Licensed under US Patent No. 7,329,412 and 11,219,706.


² A trademark of Becton, Dickinson and Company.

³ A registered trademark of CareFusion or one of its subsidiaries.

⁴ A registered trademark of 3M Company.

⁵ A registered trademark of C. R. Bard, Inc.

Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.