

ARROW ENDURANCE™ Extended Dwell Peripheral Catheter System

EXT DWELL CATH DEVICE 18GAX8CM

SKU / Article #: EDC-00818

Highlights

Sharps Safety

Contraindications & Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 20



Representative photo(s): actual product specifications may vary from that shown in the photo

ARROW ENDURANCE™ Extended Dwell Peripheral Catheter System

SKU / Article #: EDC-00818

Quantity	Component Description
1	Catheter System: 18 Ga. (1.30 mm) x 80 mm (3.15") ARROW ENDURANCE™ over 20 Ga. (1.0 mm) HW Introducer Needle with Integral 0.010" (0.26 mm) Straight Soft Tip Wire Guide
1	Label: Venous/Arterial Placement

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for

Use available under the Documents section on this page, where an IFU is available.

© 2021 Teleflex Incorporated. All Rights Reserved.