



Arrow™ Endurance™ Extended Dwell
Peripheral Catheter System
EXT DWELL CATH DEVICE 22 GA X 8 CM

SKU / Article #: EDC-00822-M


Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.
‡‡MR Safe

Sales Quantity/Case: 5







Representative photo(s): actual product specifications may vary from that shown in the photo

Arrow™ Endurance™ Extended Dwell Peripheral Catheter System SKU / Article #: EDC-00822-M	
Quantity	Component Description
1	Catheter‡‡ System: 22 Ga. (0.91 mm) x 80 mm (3.15") Arrow™ Endurance™ over 24 Ga. (0.60 mm) TW Introducer Needle with Integral 0.010" (0.26 mm) Straight Soft Tip Wire Guide
1	10 mL Luer-Lock Pre-Filled Saline Syringe
1	3 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint
1	Label: Venous/Arterial Placement
1	Tourniquet
2	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Transducer Cover with Gel and Rubber Bands

¹ A registered trademark of CareFusion or one of its subsidiaries.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.