

Representative photo(s): actual product specifications may vary from that shown in the photo

Peripheral Emergency Infusion Device - EID™

EMERGENCY INFUSION DEVICE SET: 6 FR

SKU / Article #: EI-04060

Highlights

CE Mark


Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 10

Peripheral Emergency Infusion Device - EID™ SKU / Article #: EI-04060	
Quantity	Component Description
1	Sheath: 6 Fr. x 1-1/2" (3.81 cm) Radiopaque with Dilator over 20 Ga. TW Needle
1	Spring-Wire Guide/Tube Assembly: .025" (0.64 mm) dia. x 6-7/8" (17.5 cm)
1	Scalpel: #11

Documents



Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.