

Peripheral Emergency Infusion Device - EID™

EMERGENCY INFUSION DEVICE SET: 6 FR

SKU / Article #: EI-04060

Highlights

CE Mark


Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 10

Peripheral Emergency Infusion Device - EID™ SKU / Article #: EI-04060	
Quantity	Component Description
1	Sheath: 6 Fr. x 1-1/2" (3.81 cm) Radiopaque with Dilator over 20 Ga. TW Needle
1	Spring-Wire Guide/Tube Assembly: .025" (0.64 mm) dia. x 6-7/8" (17.5 cm)
1	Scalpel: #11

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.