

Representative photo(s): actual product specifications may vary from that shown in the photo

Pediatric Jugular Puncture

PEDIATRIC JUGULAR SET: 20 GA X 12 CM

SKU / Article #: ES-04150

Highlights

CE Mark

Considerations


Not made with natural rubber latex.

Sales Quantity/Case: 10

Pediatric Jugular Puncture  
SKU / Article #: ES-04150

Quantity	Component Description
1	One-Lumen Catheter: 20 Ga. (1.2 mm OD) x 12 cm
1	Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
1	Syringe: 3 mL Luer-Slip
1	Tissue Dilator: 4 Fr. (1.35 mm) x 7.6 cm
1	Tissue Dilator: 4 Fr. (1.35 mm) x 4.1 cm
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.