



One-Lumen CVC

CVC SET: 16 GA X 8"

SKU / Article #: ES-04300

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo

One-Lumen CVC SKU / Article #: ES-04300	
Quantity	Component Description
1	One-Lumen Catheter: 16 Ga. (1.8 mm OD) x 20 cm
1	Spring-Wire Guide: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp

Documents

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.