



Central Venous Catheterization Set with Blue FlexTip[®] Catheter

CVC SET: 5 FR X 20 CM

SKU / Article #: ES-04301-M

Considerations



Not made with natural rubber latex.

Sales Quantity/Case: 10

Central Venous Catheterization Set with Blue FlexTip [®] Catheter	
SKU / Article #: ES-04301-M	
Quantity	Component Description
1	Indwelling Catheter: 5 Fr. (1.8 mm OD) x 20 cm Radiopaque Polyurethane with Blue FlexTip [®] , Integral Extension Line, Extension Line Clamps, Integral Suture Wing
1	Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 60 cm (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 18 Ga. x 6.35 to 7.20 cm RW
1	Pressure Transduction Probe
1	Arrow [®] Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator
1	Fastener: Catheter Clamp
1	Clamp: Catheter

Representative photo(s): actual product specifications may vary from that shown in the photo

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.