

## Pediatric Two-Lumen CVC

### CVC SET: 2-LUMEN 4 FR X 13 CM

SKU / Article #: ES-14402

### Highlights

CE Mark

### Considerations



Not made with natural rubber latex.

**Sales Quantity/Case: 5**

Pediatric Two-Lumen CVC SKU / Article #: ES-14402	
Quantity	Component Description
1	Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 13 cm
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 22 Ga. x 1-1/16" (2.67 cm) TW
1	Syringe: 3 mL Luer-Slip
1	Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
1	Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener

*Representative photo(s): actual product specifications may vary from that shown in the photo*

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.