



*Representative photo(s): actual product specifications may vary from that shown in the photo*

## Cath-Gard® Catheter Contamination Shield for use with 4 and 5 Fr Catheters and Arrow Percutaneous Sheath Introducer Systems

### CATH-GARD

SKU / Article #: EU-00001-ST

#### Highlights

CE Mark

#### Considerations

Not made with natural rubber latex.



Sales Quantity/Case: 10

Cath-Gard® Catheter Contamination Shield for use with 4 and 5 Fr Catheters and Arrow Percutaneous Sheath Introducer Systems

SKU / Article #: EU-00001-ST

Quantity	Component Description
1	Cath-Gard® Catheter Contamination Shield: 31-1/2" (80 cm) with Tuohy-Borst Adapter and TwistLock™ Adapter

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.