

Insertion Set for use with Arrow® 4.5 Fr. PICC

ACCESS SET

SKU / Article #: EU-00745-45

Highlights

CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5






Representative photo(s): actual product specifications may vary from that shown in the photo

Insertion Set for use with Arrow® 4.5 Fr. PICC

SKU / Article #: EU-00745-45

Quantity	Component Description
1	GlideThru™ Peel-Away Sheath: 4.5 Fr. x 2-3/4" (7 cm) Radiopaque over 4.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	SharpsAway® II Locking Disposal Cup
1	Safety Scalpel: #11

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.