



Arrow® Midline Kit

MIDLINE CATHETERIZATION KIT

SKU / Article #: EU-02041-ML

Highlights

Sharps Safety

Considerations


Not made with natural rubber latex.

Sales Quantity/Case: 5

Arrow® Midline Kit SKU / Article #: EU-02041-ML	
Quantity	Component Description
1	Midline Catheter, Trimmable: 4 Fr. (1.4 mm OD) x 20 cm, T-Port Connector, Contamination Guard and Placement Wire
1	GlideThru™ Peel-Away Sheath: 4 Fr. x 2-3/4" (7 cm) Radiopaque over 4 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	SharpsAway® II Locking Disposal Cup
1	Catheter Trimmer
1	Safety Scalpel: #11
1	Paper Tape Measure
1	Tourniquet
1	Dressing: GRIP-LOK® ¹ Catheter Fixation Device
1	Valve: MicroCLAVE™ ² Neutral Displacement Connector

¹ A registered trademark of TIDI.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.