

## One-Lumen CVC

**CVC SET: 22GA X 10CM**

SKU / Article #: EU-04522-N

### Highlights

Sharps Safety  
CE Mark

### Considerations




Not made with natural rubber latex.

**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

One-Lumen CVC	
SKU / Article #: EU-04522-N	
Quantity	Component Description
1	One-Lumen Catheter: 22 Ga. (0.9 mm OD) x 10 cm
1	Spring-Wire Guide, Nitinol, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends) with Arrow Advancer
1	Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 4 Fr. (1.35 mm) x 7.6 cm
1	Tissue Dilator: 4 Fr. (1.35 mm) x 4.1 cm
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	SharpsAway® II Locking Disposal Cup
1	Safety Scalpel: #11

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.