



## One-Lumen CVC

CVC SET: 14GA X 30CM

SKU / Article #: EU-04730-CVT

### Highlights

CE Mark

### Considerations

Not made with natural rubber latex.

**Sales Quantity/Case: 10**

One-Lumen CVC SKU / Article #: EU-04730-CVT	
Quantity	Component Description
1	One-Lumen Catheter: 14 Ga. (2.2 mm OD) x 30 cm, Side Holes
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
1	ECG Extension Cable
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Scalpel: #11

*Representative photo(s): actual product specifications may vary from that shown in the photo*

