

## Pressure Injectable One-Lumen PICC

PI PICC KIT: 4 FR X 55 CM

SKU / Article #: EU-05541-HPMSB

### Highlights

Sharps Safety  
CE Mark

### Considerations

Not made with natural rubber latex.

**Sales Quantity/Case: 5**

Pressure Injectable One-Lumen PICC SKU / Article #: EU-05541-HPMSB	
Quantity	Component Description
1	One Lumen Catheter: 4 Fr. (1.45 mm OD) x 55 cm, Pressure Injectable, T-Port Connector, Blue FlexTip®, Contamination Guard and Placement Wire
1	GlideThru™ Peel-Away Sheath: 4 Fr. x 2-3/4" (7 cm) Radiopaque over 4 Fr. Dilator
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
1	Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
1	Injection Needle: Eclipse™ <sup>1</sup> 23 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Injection Needle: Eclipse™ <sup>1</sup> 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Syringe: 10 mL Luer-Lock
1	Syringe: 10 mL Luer-Slip
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp

*Representative photo(s): actual product specifications may vary from that shown in the photo*

## Pressure Injectable One-Lumen PICC

SKU / Article #: EU-05541-HPMSB





Quantity	Component Description
1	SecondSite™ Adjustable Hub: Fastener
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations
1	Towel
1	Safety Scalpel: #11
1	Patient ID Card
1	Sterile Procedure Sign
1	Patient Information Booklet
1	Chart Sticker
1	PICC Line Insertion Checklist
2	Paper Tape Measure
1	Tourniquet
5	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask with Eye Shield
1	Dressing: STATLOCK® <sup>2</sup> PICC Plus Catheter Stabilization Device
1	Surgical Apparel: Mask
1	Tape: Steri-Strip® <sup>3</sup>
1	HemoHopper® Fluid Receptacle

<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of C. R. Bard, Inc.

<sup>3</sup> A registered trademark of 3M Company.

# Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.