

Representative photo(s): actual product specifications may vary from that shown in the photo

Three-Lumen CVC for Femoral Vein Insertion

CVC SET: 3-LUMEN 5.5FR X 30CM

SKU / Article #: EU-14553-N

Highlights

Sharps Safety  
CE Mark


Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

Three-Lumen CVC for Femoral Vein Insertion SKU / Article #: EU-14553-N	
Quantity	Component Description
1	Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 30 cm
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
1	Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
1	Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	Safety Scalpel: #11

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.