

NO IMAGE
AVAILABLE

Pressure Injectable Three-Lumen CVC
PI CVC KIT: 3-LUMEN 7 FR X 20 CM

SKU / Article #: EU-15703-HPSKM

Highlights

Sharps Safety
CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5


Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Three-Lumen CVC SKU / Article #: EU-15703-HPSKM	
Quantity	Component Description
1	Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm, Pressure Injectable
1	Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer with ECG Mark
1	Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Injection Needle: Eclipse™1 25 Ga. x 1" (2.54 cm)
1	Pressure Transduction Probe
1	Syringe: 3 mL Luer-Lock
1	Syringe: 10 mL Luer-Slip
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 8 Fr. (2.7 mm) x 9 cm
3	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener

Pressure Injectable Three-Lumen CVC SKU / Article #: EU-15703-HPSKM	
Quantity	Component Description
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	Towel: 13" x 18" (33 cm x 46 cm)
1	Maximal Barrier Drape: 55" x 94-1/2" (140 cm x 240 cm) with 4" (10 cm) fenestration, with adhesive and covered perforation
1	Safety Scalpel: #11
3	Prep Sponge Swab
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm® ² I.V. Advanced 3-1/4" x 4-1/2" (8.5 cm x 11.5 cm)
1	Dressing: GRIP-LOK® ³ Catheter Fixation Device
1	Surgical Apparel: Cap
1	Surgical Apparel: Impervious Gown
1	Surgical Apparel: Mask with Shield
1	Transducer Cover: 4-12/16" x 95" (12 cm x 244 cm)
1	Tray: Prep

¹ A trademark of Becton, Dickinson and Company.
² A registered trademark of 3M Company.
³ A registered trademark of TIDI.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for

Use available under the Documents section on this page, where an IFU is available.