

Five-Lumen CVC

CVC SET: 5-LUMEN 9.5FR X 20CM

SKU / Article #: EU-15955-CVET

Highlights

CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5



Representative photo(s): actual product specifications may vary from that shown in the photo

Five-Lumen CVC	
SKU / Article #: EU-15955-CVET	
Quantity	Component Description
1	Five-Lumen Catheter: 9.5 Fr. (3.3 mm OD) x 20 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
1	ECG Extension Cable
5	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Scalpel: #11

