

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Arrowg+ard Blue®
Two-Lumen CVC

CVC SET: 2-LUMEN 7FR X 16CM

SKU / Article #: EU-22702-HPS

Highlights

- Antimicrobial: Arrowg+ard Blue*
- Sharps Safety
- CE Mark




Considerations

- *Chlorhexidine
- Not made with natural rubber latex.

Sales Quantity/Case: 10

Pressure Injectable Arrowg+ard Blue® Two-Lumen CVC SKU / Article #: EU-22702-HPS	
Quantity	Component Description
1	Arrowg+ard Blue® Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm, Pressure Injectable
1	Spring-Wire Guide, Nitinol, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 8 Fr. (2.7 mm) x 9 cm
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	Safety Scalpel: #11

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.