




## Pressure Injectable One-Lumen PICC

### PICC KIT: 4 FR X 50 CM

SKU / Article #: EU-25041-HPMSB

### Highlights

Sharps Safety  
CE Mark

### Considerations

Not made with natural rubber latex.

### Sales Quantity/Case: 5

*Representative photo(s): actual product specifications may vary from that shown in the photo*

| Pressure Injectable One-Lumen PICC<br>SKU / Article #: EU-25041-HPMSB |                                                                                                                                                             |
|-----------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Quantity                                                              | Component Description                                                                                                                                       |
| 1                                                                     | One Lumen Catheter: 4 Fr. (1.45 mm OD) x 50 cm, Pressure Injectable, T-Port Connector, Blue FlexTip®, Contamination Guard and Placement Wire                |
| 1                                                                     | GlideThru™ Peel-Away Sheath: 4 Fr. x 2-3/4" (7 cm) Radiopaque over 4 Fr. Dilator                                                                            |
| 1                                                                     | Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer |
| 1                                                                     | Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip                                                                                             |
| 1                                                                     | Injection Needle: Eclipse™ <sup>1</sup> 23 Ga. x 1" (2.54 cm)                                                                                               |
| 1                                                                     | Injection Needle: Eclipse™ <sup>1</sup> 25 Ga. x 1" (2.54 cm)                                                                                               |
| 1                                                                     | Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW                                                                                                      |
| 2                                                                     | Syringe: 10 mL Luer-Lock                                                                                                                                    |
| 1                                                                     | Syringe: 3 mL Luer-Lock                                                                                                                                     |
| 1                                                                     | SecondSite™ Adjustable Hub: Catheter Clamp                                                                                                                  |
| 1                                                                     | SecondSite™ Adjustable Hub: Fastener                                                                                                                        |
| 1                                                                     | SharpsAway® II Locking Disposal Cup                                                                                                                         |

# Pressure Injectable One-Lumen PICC

SKU / Article #: EU-25041-HPMSB


| Quantity | Component Description                                                             |
|----------|-----------------------------------------------------------------------------------|
| 1        | SharpsAway® Disposal Cup                                                          |
| 1        | Catheter Trimmer                                                                  |
| 1        | Towel: 24" x 36" (61 cm x 91 cm)                                                  |
| 1        | Drape: 60" x 76" (152 cm x 193 cm)                                                |
| 1        | Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration        |
| 1        | Towel                                                                             |
| 1        | Safety Scalpel: #11                                                               |
| 1        | Patient ID Card                                                                   |
| 1        | Sterile Procedure Sign                                                            |
| 1        | Patient Information Booklet                                                       |
| 1        | Chart Sticker                                                                     |
| 1        | PICC Line Insertion Checklist                                                     |
| 1        | Paper Tape Measure                                                                |
| 1        | Tourniquet                                                                        |
| 5        | Gauze Pad: 2" x 2" (5 cm x 5 cm)                                                  |
| 10       | Gauze Pad: 4" x 4" (10 cm x 10 cm)                                                |
| 1        | Dressing: Tegaderm® <sup>2</sup> I.V. Advanced 3-1/4" x 4-1/2" (8.5 cm x 11.5 cm) |
| 1        | Dressing: GRIP-LOK® <sup>3</sup> Catheter Fixation Device                         |
| 1        | Surgical Apparel: Cap                                                             |
| 1        | Surgical Apparel: Impervious Gown                                                 |
| 1        | Dressing: Steri-Strip® <sup>2</sup> 1/2" x 4" (1.2 cm x 10 cm)                    |
| 1        | Surgical Apparel: Mask                                                            |
| 1        | Surgical Apparel: Mask with Shield                                                |
| 1        | Transducer Cover: 4-12/16" x 95" (12 cm x 244 cm)                                 |
| 1        | Valve: MicroCLAVE™ <sup>4</sup> Neutral Displacement Connector                    |
| 1        | HemoHopper® Fluid Receptacle                                                      |

<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of 3M Company.

<sup>3</sup> A registered trademark of TIDI.  
<sup>4</sup> A trademark of ICU Medical, Inc.

Documents



SCROLL TO VIEW MORE ➔

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.