

Pressure Injectable One-Lumen PICC

PICC KIT: 4 FR X 55 CM

SKU / Article #: EU-25541-IR

Highlights

Sharps Safety
CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

NO IMAGE
AVAILABLE

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable One-Lumen PICC





SKU / Article #: EU-25541-IR

Quantity	Component Description
1	One Lumen Catheter: 4 Fr. (1.45 mm OD) x 55 cm, Pressure Injectable, Blue FlexTip®
1	GlideThru™ Peel-Away Sheath: 4 Fr. x 2-3/4" (7 cm) Radiopaque over 4 Fr. Dilator
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
1	Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 51-3/16" (130 cm)(Straight Soft Tip on One End - Straight Stiff Tip on Other)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	SharpsAway® II Locking Disposal Cup
1	Catheter Trimmer
1	Safety Scalpel: #11
1	Patient ID Card
1	Patient Information Booklet
1	Chart Sticker
1	Paper Tape Measure
1	Tourniquet
1	Dressing: GRIP-LOK® ¹ Catheter Fixation Device
1	Valve: MicroCLAVE™ ² Neutral Displacement Connector

¹ A registered trademark of TIDI.

² A trademark of ICU Medical, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.