

## Arrowgard Blue® Three-Lumen CVC

### CVC SET: 3-LUMEN 5.5FR X 8CM

SKU / Article #: EU-25553-N


### Highlights

Antimicrobial: Arrowgard Blue\*  
 Sharps Safety  
 CE Mark

### Considerations

\*Chlorhexidine  
 Not made with natural rubber latex.




**Sales Quantity/Case: 5**



Representative photo(s): actual product specifications may vary from that shown in the photo

Arrowgard Blue® Three-Lumen CVC SKU / Article #: EU-25553-N	
Quantity	Component Description
1	Arrowgard Blue® Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 8 cm
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
1	Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
1	Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	Safety Scalpel: #11

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.