

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Three-Lumen PICC

PICC KIT: 6 FR X 55 CM

SKU / Article #: EU-25563-HP

Highlights

Sharps Safety
CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

Pressure Injectable Three-Lumen PICC





SKU / Article #: EU-25563-HP

| Quantity | Component Description |
|----------|---|
| 1 | Three Lumen Catheter: 6 Fr. (2.14 mm OD) x 55 cm, Pressure Injectable, T-Port Connector, Blue FlexTip®, Contamination Guard and Placement Wire |
| 1 | GlideThru™ Peel-Away Sheath: 6 Fr. x 2-3/4" (7 cm) Radiopaque over 6 Fr. Dilator |
| 1 | Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer |
| 1 | Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW |
| 1 | Syringe: 10 mL Luer-Lock |
| 1 | SecondSite™ Adjustable Hub: Fastener |
| 1 | SecondSite™ Adjustable Hub: Catheter Clamp |
| 1 | SharpsAway® II Locking Disposal Cup |
| 1 | Catheter Trimmer |
| 1 | Safety Scalpel: #11 |
| 1 | Patient ID Card |
| 1 | Patient Information Booklet |
| 1 | Chart Sticker |
| 1 | Paper Tape Measure |
| 1 | Tourniquet |
| 1 | Dressing: GRIP-LOK® ¹ Catheter Fixation Device |
| 3 | Valve: MicroCLAVE™ ² Neutral Displacement Connector |

¹ A registered trademark of TIDI.

² A trademark of ICU Medical, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

© 2025 Teleflex Incorporated. All Rights Reserved.