

Arrowgard Blue® Three-Lumen CVC

CVC SET: 3-LUMEN 5.5FR X 13CM

SKU / Article #: EU-26553-N

Highlights

- Antimicrobial: Arrowgard Blue*
- Sharps Safety
- CE Mark

Considerations




- *Chlorhexidine
- Not made with natural rubber latex.

Sales Quantity/Case: 5

Arrowgard Blue® Three-Lumen CVC	
SKU / Article #: EU-26553-N	
Quantity	Component Description
1	Arrowgard Blue® Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 13 cm
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
1	Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
1	Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	Safety Scalpel: #11

Representative photo(s): actual product specifications may vary from that shown in the photo

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.