

Representative photo(s): actual product specifications may vary from that shown in the photo

Arrowg+ard Blue® Two-Lumen CVC

CVC SET: 2-LUMEN 7FR X 20CM

SKU / Article #: EU-26702-CVT

Highlights

Antimicrobial: Arrowg+ard Blue*
CE Mark


Considerations

*Chlorhexidine
Not made with natural rubber latex.

Sales Quantity/Case: 10

Arrowg+ard Blue® Two-Lumen CVC SKU / Article #: EU-26702-CVT	
Quantity	Component Description
1	Arrowg+ard Blue® Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
1	ECG Extension Cable
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Scalpel: #11

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.