



Pressure Injectable Two-Lumen Midline

PI MIDLINE 2-L: 5FR X 20CM

SKU / Article #: EU-32052

Highlights

Sharps Safety
CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5



Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Two-Lumen Midline SKU / Article #: EU-32052	
Quantity	Component Description
1	Two-Lumen TaperFree® Catheter†††: 5 Fr.(1.81 mm OD) x 20 cm Pressure Injectable Radiopaque Polyurethane with T-Port Connector, Extension Line Clamps, Contamination Guard and Stiffening Stylet
1	GlideThru™ Peel-Away Sheath: 5 Fr. x 2-3/4" (7 cm) Radiopaque over 5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	SharpsAway® II Locking Disposal Cup
1	Catheter Trimmer
1	Safety Scalpel: #11
1	Paper Tape Measure
1	Tourniquet
1	Dressing: GRIP-LOK® ¹ Catheter Fixation Device
2	Valve: MicroCLAVE™ ² Neutral Displacement Connector

¹ A registered trademark of TIDI.

² A trademark of ICU Medical, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.