

Arrowgard Blue Plus® Three-Lumen CVC

CVC KIT: 3-LUMEN 7FR X 20CM

SKU / Article #: EU-45703-MSB

Highlights

Antimicrobial: Arrowgard Blue Plus*
 Sharps Safety
 CE Mark

Considerations

*Chlorhexidine
 Not made with natural rubber latex.

Sales Quantity/Case: 5

| Arrowgard Blue Plus® Three-Lumen CVC | |
|--------------------------------------|---|
| SKU / Article #: EU-45703-MSB | |
| Quantity | Component Description |
| 1 | Arrowgard Blue Plus® ¹ Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm |
| 1 | Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark |
| 1 | Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle |
| 1 | Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip |
| 1 | Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW |
| 1 | Injection Needle: Eclipse ^{TM2} 23 Ga. x 1" (2.54 cm) |
| 1 | Injection Needle: Eclipse ^{TM2} 25 Ga. x 1" (2.54 cm) |
| 1 | Pressure Transduction Probe |
| 1 | Syringe: 5 mL Luer-Lock |
| 1 | Syringe: 3 mL Luer-Lock |
| 1 | Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL |

Representative photo(s): actual product specifications may vary from that shown in the photo

Arrowg+ard Blue Plus Three-Lumen CVC




SKU / Article #: EU-45703-MSB

| Quantity | Component Description |
|----------|--|
| 1 | Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm |
| 1 | SecondSite™ Adjustable Hub: Fastener |
| 1 | SecondSite™ Adjustable Hub: Catheter Clamp |
| 1 | SharpsAway® II Locking Disposal Cup |
| 1 | SharpsAway® Disposal Cup |
| 1 | Towel: 13" x 18" (33 cm x 46 cm) |
| 1 | Maximal Barrier Drape: 55" x 94-1/2" (140 cm x 240 cm) with 4" (10 cm) fenestration, with adhesive and covered perforation |
| 1 | Needle Holder |
| 1 | Safety Scalpel: #11 |
| 3 | Stopcock: 4-way High-Flow |
| 1 | Sterile Procedure Sign |
| 2 | Prep Sponge Swab |
| 2 | Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 6 | Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1 | Surgical Apparel: Cap |
| 1 | Surgical Apparel: Impervious Gown |
| 1 | Surgical Apparel: Mask with Shield |
| 1 | Suture: (0) Silk with Curved Needle |
| 1 | Tray: Prep |

¹ Licensed under US Patent No. 7,329,412.

² A trademark of Becton, Dickinson and Company.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.