

## Pressure Injectable Arrowg+ard Blue Plus® Two-Lumen CVC

PI CVC SET: 2-LUMEN 8FR X 20CM

SKU / Article #: EU-45802-HPS

### Highlights

Antimicrobial: Arrowg+ard Blue Plus\*  
 Sharps Safety  
 CE Mark

### Considerations

\*Chlorhexidine  
 Not made with natural rubber latex.




**Sales Quantity/Case: 10**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Pressure Injectable Arrowg+ard Blue Plus® Two-Lumen CVC	
SKU / Article #: EU-45802-HPS	
Quantity	Component Description
1	Arrowg+ard Blue Plus® <sup>1</sup> Two-Lumen Catheter: 8 Fr. (2.8 mm OD) x 20 cm, Pressure Injectable
1	Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 9 Fr. x 10.2 cm
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	Safety Scalpel: #11

<sup>1</sup> Licensed under US Patent No. 7,329,412.

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

© 2026 Teleflex Incorporated. All Rights Reserved.